

HOT PLATES

HOT PLATE

30 x 30 cm

M 3030

Temperature working range	50 °C to 300 °C
Temperature control	Digital thermostat
Thermostat setting increments	in 1 °C steps
Thermostat working accuracy	± 1 °C
Power rating	1500 W
Fuse	25 Amps.
Casing material	Sheet aluminum finished in stoved teflon paint
Heating plate material	Teflon coated sheet aluminum
Heating plate dimensions	30 x 30 cm
Power supply	AC 190 - 240 V / 50-60 Hz
Overall dimensions (w x h x d)	30 x 14 x 30 cm
Net weight	8 kg
Gross weight	9 kg

